

Rotograph EVO 3D

DAP Dose per Area Product					
panoramic exam			unit of measure: $\mu\text{Gy}\cdot\text{m}^2$		
mA	6	7	8	9	10
kV					
60	5.0	5.8	6.6	7.4	8.2
62	5.2	6.1	7.0	7.9	8.8
64	5.6	6.6	7.5	8.4	9.4
66	6.1	7.1	8.0	9.0	10.0
68	6.4	7.5	8.6	9.6	10.7
70	6.9	8.1	9.2	10.3	11.5
72	7.2	8.4	9.6	10.8	11.9
74	7.6	8.8	10.1	11.4	12.6
76	8.0	9.3	10.6	11.9	13.3
78	8.3	9.8	11.1	12.5	14.0
80	8.8	10.3	11.8	13.3	14.7
82	9.1	10.6	12.1	13.6	15.1
84	9.6	11.2	12.8	14.4	16.0
86	10.1	11.8	13.5	15.1	16.8